

Datasheet

Catalog Number: NBL-243150

Version: RN6.0



Product Name	NebuSelect™ Anti-mGluR1 Polyclonal Antibody
Alias/Synonyms	GRM1, GRM1A, GPRC1A, MGlu1, MGLUR1, SCAR13, MGLUR1A
Host / Isotype	Rabbit / IgG
Type	Polyclonal
Appearance	Liquid
Description	NebuSelect™ Anti-mGluR1 Polyclonal antibody (NBL-243150) is a highly specific antibody directed against the GluR1 protein. The antibody can be used in Immunofluorescence and immunocytochemistry applications, and recognizes mGluR1 in human and mouse samples.
Tested Applications	IHC, IF-P, ELISA
Tested Reactivity	human, mouse
Immunogen	A synthetic peptide corresponding to the center region of the Human GluR1 protein.
Buffer	0.2 µm filtered solution in PBS
Size	100ul
Storage	This product remains stable up to 18 months at -20°C. Avoid repeated freeze-thaw cycles.