

Datasheet

Catalog Number: NBL-243151

Version: RN6.0



Product Name	NebuSelect™ Anti-mGluR5 Polyclonal antibody
Alias/Synonyms	GRM5, GPRC1E, MGLUR5, Mglur5b, MGlu5
Host / Isotype	Rabbit / IgG
Type	Polyclonal
Appearance	Liquid
Description	NebuSelect™ Anti-mGluR5 Polyclonal antibody (NBL-243151) is a highly specific antibody directed against mGluR5. The antibody can be used in IHC, IF and ELISA applications. It has been designed to recognize mGluR5 from rat, mouse, and human samples.
Tested Applications	IHC, IF, ELISA
Tested Reactivity	Rat (predicted in Human, Mouse)
Immunogen	Synthetic peptide derived from human GRM5: 201-300.
Buffer	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Size	100ul
Storage	This product remains stable up to 18 months at -20°C. Avoid repeated freeze-thaw cycles.