

Datasheet

Catalog Number: NBL-243156

Version: RN6.0



Product Name	NebuSelect™ Anti-CD207 Polyclonal antibody
Alias/Synonyms	Langerin, C-type lectin domain family 4 member K, CLEC4K, CD207 molecule, langerin
Host / Isotype	Rabbit / IgG
Type	Polyclonal
Appearance	Liquid
Description	NebuSelect™ Anti-CD207 Polyclonal antibody(NBL-243156) is a highly specific antibody directed against the CD207 protein. The antibody can be used in western blot and immunocytochemistry applications, and recognizes CD207 in human and mouse samples.
Tested Applications	WB, IHC WB: 1:500-1:1000 IHC: 1:1000-1:4000
Tested Reactivity	human, mouse
Immunogen	CD207 Fusion Protein
Buffer	0.2 µm filtered solution in PBS
Size	100ul; 500ul; 1ml
Storage	This product remains stable up to 18 months at -20°C. Avoid repeated freeze-thaw cycles.