

Datasheet

Catalog Number: NBL-243154

Version: RN6.0



Product Name	NebuSelect™ Anti-CD1a Polyclonal antibody
Alias/Synonyms	CD1, CD1A, CD1a molecule, FCB6, HTA1, hTa1 thymocyte antigen, R4, T6
Host / Isotype	Rabbit / IgG
Type	Polyclonal
Appearance	Liquid
Description	NebuSelect™ Anti-CD1a Polyclonal antibody(NBL-243154) is a highly specific antibody directed against the CD1a protein. The antibody can be used in immunocytochemistry applications, and recognizes CD1a in human, mouse and pig samples.
Tested Applications	IHC IHC: 1:1000-1:4000
Tested Reactivity	human, mouse, pig
Immunogen	CD1a Fusion Protein
Buffer	0.2 µm filtered solution in PBS
Size	100ul; 500ul; 1ml
Storage	This product remains stable up to 18 months at -20°C. Avoid repeated freeze-thaw cycles.