

Datasheet

Catalog Number: NBL-243155

Version: RN6.0



Product Name	NebuSelect™ Recombinant human PDGF-AA
Alias/Synonyms	Platelet-derived growth factor subunit A;PDGF subunit A;PDGF-1;Platelet-derived growth factor A chain;Platelet-derived growth factor alpha polypeptide; PDGFA;PDGF1
Appearance	Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.
Description	NebuSelect™ Recombinant human PDGF-AA(Cat#NBL-243155) is expressed in E.coli.
Accession	P04085
MW	The protein has a predicted MW of 14.1 kDa.
Purity	> 90%
Endotoxin	< 1EU/ug(LAL Method)
Size	100ug; 500ug
Storage	Short term 4°C; Long term -20°C. Avoid repeated freeze-thaw cycles.