

Datasheet

Catalog Number: NBL-243157

Version: RN6.0



Product Name	NebuSelect™ Anti-Connexin32 Monoclonal antibody
Alias/Synonyms	CX32, GJB1, Connexin32
Host / Isotype	Mouse / IgG1κ
Type	Monoclonal
Appearance	Liquid(0.5mg/mL)
Description	NebuSelect™ Anti-Connexin32 Monoclonal antibody(CAT#NBL-243157) recognizes human, mouse, rat, and chicken Connexin 32. This antibody has been reported for use in western blotting, immunohistochemical staining of frozen tissue sections, immunohistochemical staining of formalin-fixed paraffin embedded tissue sections, microscopy, and immunocytochemistry.
Tested Applications	WB, IHC(5~10μg/mL), ICC/IF
Tested Reactivity	Chicken, Human, Mouse, Rat
Immunogen	Recombinant Connexin32 protein.
Buffer	PBS, pH 7.2
Size	100ul
Storage	4°C.