

Datasheet

Catalog Number: NBL-243158

Version: RN6.0



Product Name	NebuSelect™ Anti-MTAP Monoclonal antibody
Alias/Synonyms	c86fus, MSAP, MTA phosphorylase, MTAP, MTAPase
Host / Isotype	Mouse / IgG1k
Type	Monoclonal
Dilution	WB : 1:2000-1:10000 IHC : 1:50-1:1000 IF/ICC : 1:400-1:1600
Appearance	Liquid
Description	NebuSelect™ Anti-MTAP Monoclonal antibody(CAT#NBL-243158) recognizes human, mouse, and rat MTAP. This antibody has been reported for use in western blotting, immunohistochemical staining of frozen tissue sections, immunohistochemical staining of formalin-fixed paraffin embedded tissue sections, and Immunofluorescence.
Tested Applications	WB, IHC, IF/ICC, ELISA
Tested Reactivity	human, mouse, rat
Immunogen	recombinant MTAP (BC018625, 1 a.a. ~ 154 a.a) full-length protein
Buffer	0.2 µm filtered solution in PBS, pH 7.2
Size	100ul; 500ul; 1ml
Storage	This product remains stable up to 18 months at -20°C. Avoid repeated freeze-thaw cycles.