


Datasheet

Catalog Number: NBL-250783

Version: RN6.0



Product Name	NebuSelect™ Recombinant Mouse IL-22R Protein, His tag
Alias/Synonyms	IL22RA1; CRF2-9; ZcytoR11; Interleukin 22 Receptor,Alpha 1; Cytokine receptor class-II member 9
Appearance	Freeze-dried powder from 0.22 µm filtered solution in PBS, pH7.4 with 5% Trehalose and 5% Proclin300.
Description	NebuSelect™ Recombinant Mouse IL-22R Protein, His tag (Cat#NBL-250783) is expressed in E.coli.
AA Sequence	Pro258-Cys552
MW	The protein has a predicted MW of 40 kDa.
Activity	Not validated for activity
Purity	> 90%
Endotoxin	< 0.1EU/ug(LAL Method)
Size	100ug; 500ug; 1mg
Storage	Short term 4°C; Long term -20°C. Avoid repeated freeze-thaw cycles.
Application	Positive Control; Immunogen; SDS-PAGE; WB.
Note	For research use only.



M: Marker; Lane 1: NBL250783

4-20% SDS-PAGE Coomassie Staining (2ug)