

Datasheet

Catalog Number: NBL-250780

Version: RN6.0



Product Name	NebuSelect™ Recombinant Human HSA
Alias/Synonyms	Serum Albumin; HAS; ALB
Appearance	Liquid from 0.22 µm filtered solution in PBS, pH7.4 with 5% Trehalose and 5% Mannitol.
Description	NebuSelect™ Recombinant Human HSA(Cat#NBL-246739) is expressed in HEK 293 cells.
Accession	P02768
AA Sequence	A25-L609
MW	The protein has a predicted MW of 66.5 kDa.
Activity	Not validated for activity
Purity	> 95%
Endotoxin	< 0.1EU/ug(LAL Method)
Size	100ul; 500ul; 1ml
Storage	Short term 4°C; Long term -20°C. Avoid repeated freeze-thaw cycles.
Application	Positive Control; Immunogen; SDS-PAGE; WB.
Note	For research use only.